




The ERGYS Data Management System

# Physician's Electronic Prescription Station



Prescription date

Patient name

Physician name

Days per week for ERGYS therapy (1-7)

Max session time (3 to 60 minutes)

Max run time (3 to 60 minutes)

Rest time (3 to 60 minutes)

Workload (in 1/8 Kp increments)

manual resistance

automatic resistance

1:ERGYS1+/-450	2:SINE 300	3:SINE 400	4:SINE 500	5:CAMELBACK 500
6:SINE-50- 300	7:SINE-100- 350	8:SINE-100- 400	15:ERGYS1 450	

Pulses per second (30, 40, 50, 60)

Control RPM Setting:

10 14 18 21 25 28 32 35 39 42 46 49 53 56 60 63 67 70

Fatigue RPM Setting:

10 14 18 21 25 28 32 35 39 42 46 49 53 56 60 63 67 70

Max Stim settings (0 to 140 MA):

RQUAD          LQUAD          RGLUT          LGLUT          RHAMS          LHAMS

Allow ERGYS rider to change Waveshape, RPM, and Max Stim Parameters:    YES    NO

Therapeutic Alliances Inc.  
333 North Broad Street  
Fairborn, Ohio 45324 USA



(937) 879-0734 • (937) 879-5211 (fax) • www.ERGYS.com • info@ERGYS.com

ERGYS® is a trademark of Therapeutic Alliances Inc.